Title: Squash Preparation of Grasshopper Testis and Study of the Various Stages of Meiosis

Objective: To prepare a temporary stained squash of grasshopper testis and observe the various stages of meiosis under a microscope.

Introduction: Meiosis is a specialized type of cell division that reduces the chromosome number by half, producing haploid gametes. It consists of two successive divisions: **Meiosis I** (reductional division) and **Meiosis II** (equational division). The grasshopper testis is an ideal material for studying meiosis due to its large, easily distinguishable chromosomes.

Materials Required:

- 1. Live grasshopper (freshly dissected testis)
- 2. Saline solution (0.9% NaCl)
- 3. Hydrochloric acid (HCl, 1N)
- 4. Aceto-orcein or Acetocarmine stain
- 5. Distilled water
- 6. Watch glass
- 7. Scalpel or fine scissors
- 8. Forceps
- 9. Glass slides and cover slips
- 10. Filter paper or blotting paper
- 11. Compound light microscope

Procedure:

1. Dissection and Collection of Testis:

- o Anesthetize the grasshopper and carefully dissect it to expose the testis.
- Remove the testis using fine forceps and place it in a watch glass containing saline solution.

2. Fixation and Hydrolysis:

- Transfer the testis to a watch glass containing 1N HCl and incubate for 3–5 minutes to soften the tissue.
- o Rinse the tissue with distilled water and blot dry with filter paper.

3. Staining:

- o Transfer the testis onto a clean glass slide.
- Add a few drops of aceto-orcein or acetocarmine stain and allow it to stain for 5–10 minutes.

4. Squashing:

- o Cover the stained testis with a cover slip.
- o Gently press the cover slip with the blunt end of a forceps or the back of a pencil to spread the cells into a thin layer.
- o Remove excess stain using filter paper.

5. Microscopic Observation:

- o Place the prepared slide on the microscope stage.
- o Start with a low-power objective (10x) to locate the dividing cells.
- Switch to high power (40x) to observe the different stages of meiosis.

Observations:

• Meiosis I:

- o *Prophase I:* Chromosomes condense, homologous chromosomes pair (synapsis), and crossing over occurs.
- o *Metaphase I:* Homologous chromosome pairs align at the equatorial plate.
- Anaphase I: Homologous chromosomes separate and move to opposite poles.
- o Telophase I: Two haploid daughter cells are formed.

• Meiosis II:

- o *Prophase II:* Chromosomes condense again without DNA replication.
- o Metaphase II: Chromosomes align at the equatorial plate.
- o Anaphase II: Sister chromatids separate and move to opposite poles.
- o Telophase II: Four haploid gametes are formed.

Discussion: The experiment successfully demonstrated the different stages of meiosis in grasshopper testis cells. Hydrolysis with HCl softened the tissues, and staining enhanced chromosome visibility. Meiosis I reduced chromosome number, ensuring genetic diversity, while Meiosis II resulted in haploid gametes. The observed stages confirm the importance of meiosis in gametogenesis and genetic variation.

Conclusion: The squash preparation technique enabled clear observation of meiotic stages in grasshopper testis cells. This experiment highlights the fundamental role of meiosis in sexual reproduction and chromosome reduction.




Figure 1: Different stages of meiosis in grasshopper testis squash